


BULLETIN

NO.: 20-
DATE: May 15, 2020
TO: All Personnel
FROM: Chris Heiser, Deputy Chief, Emergency Medical Services
SUBJECT: Gibeck Humid Vent Filter Compact

Effective immediately, personnel shall discontinue using the Gibeck Humid Vent Compact filter with any in-line nebulizer.

The manufacturer has recommended to discontinue using of the Humid Vent filter during any in-line nebulization because it may degrade the filter's performance and require replacement.



The Humid Vent filter is still effective when using a BVM without an in-line nebulizer.



Contact your chain of command with questions or concerns.